




Wearable Wireless Breast Pump

User Manual



Read this user manual carefully before using your electric breast pump, and save the user manual for future reference.


Overview



Parts list

- ① Silicone shield x1
- ② Pump motor x1
- ③ Membrane x1
- ④ Valve x1
- ⑤ 150ml Cup x1
- ⑥ Dust cover x1
- ⑦ Silicone cylinder x1

Control panel



- ① Display
- ② Change Mode
- ③ Decrease
- ④ Increase
- ⑤ ON / OFF
- ⑥ Interface

Specification


Rated Voltage	5V
Rated Power	3W
Ampere:	1A
Battery Capacity:	900mAh
Power Adapter:	5V- 1A
Auto Power Off:	30 minutes

Assembling

- ① Install the valve into the bottom of the connector, making sure it is in the correct orientation.
- ② Place the assembled connector into the milk collector, aligning the "ear" part of the connector with the internal groove of the milk collector, and press the disc part of the connector with both hands to secure the edges firmly to ensure a complete seal.
- ③ Install the silicone shield on the milk collector and align the guide plate on the silicone shield with the milk collector. Align the deflector on the top and press the edges with both hands to secure it in place.



- ④ Press the pump motor onto the milk collector, making sure it fits snugly.
- ⑤ Place the silicone shield into the milk collector and press it with your hands to tighten the edges to ensure a good seal.
- ⑥ Remove the dust cover and connect the USB cable to charge the pump. The charging process takes 130 minutes and the pump can be used for about 90 minutes.




Note

The power adapter is a common standard accessory on the market. This product is not provided. Power adapter specification: 5V = 1A.

TO begin pumping

- ① Press the silicone shield firmly against the breast, leave no gaps and make sure your nipples are centered. Keep it upright. Be careful not to distort the position of the product and the chest, otherwise the breast pump will not have suction.



- ② Lightly press the power button to start the breast pump, and the timing display turns on.



- ③ Lightly press the "M" key to switch the mode, the display shows "9" for the breast pumping mode. There are 9 levels in this mode. You can select the suction power that makes you comfortable by touching the "+" and "-" keys. During postpartum breastfeeding, breast discomfort such as swelling and lumps are prone to occur. Using this mode can alleviate breast problems.



- ④ Press the "M" key again to switch to the massage mode, the display shows "9" for the massage mode. There are 9 levels in this mode. Each level have 9 short suction and short release plus 1 long suction and long release as a cycle. You can press the "+" and "-" keys to select the suction power that makes you comfortable. In the early postpartum period, the breast duct is not easy to conduct, and the breast can be stimulated by changing the mode to promote breast milk secretion.



Tips: The machine has a memory function. When you turn it on again, the machine will automatically enter the working mode and suction position you used last time.

- ⑤ When you have finished pumping, press the power button to turn off the breast pump, then carefully remove the breast pump from breast.



Cleaning and sterilizing

- For the first time and before each use, you must disassemble the pump to wash and sanitize all the parts that come into contact with the breast milk.
- Never wash or sanitize the pump motor, You can wipe this part with a clean, soft cloth.
- Clean the product with a bottle brush using warm water and mild detergent.
- It can be sanitized by putting it in hot water for 5 minutes, but be careful not to add the parts to the water while it is boiling. Let the water calm down first, then add the parts to sanitize.
- Take care when cleaning the valve and silicone diaphragm, as damage to these parts will prevent the breast pump from functioning properly.
- To clean the valve, rub it gently between your fingers in warm, soapy water. Do not insert any objects as this may cause damage.

Sitting position for suction


Relax the body, sit in a comfortable chair and lean slightly forward (use cushion to support back). Make sure to keep the breast pump and milk collector upright.



Bra adjustment buckle


Adjust the length of bra straps, make suction more comfortable.

One end of the bra adjustment buckle hung on the inside buckle of nursing bra. The length of the bra strap can be adjusted by hanging the nursing bra buckle on the other end of the bra adjustment buckle.



Pouring milk

Please keep breast pump upright, and remove the pump motor from pump body carefully, to pour the milk into the bottle, make sure the gap of milk collector against the bottle mouth.



Storing your breast milk

Storing breast milk

Store the breast milk in clean and sterilized storage bag or sealed bottle, and put them in fridge immediately.

Unfreeze breast milk

Put the frozen breast milk in the fridge to unfreeze overnight. In case of any emergency, you can defrost the milk in a bowl of hot water, make sure the bottle stand up, or bottle warmer to heat the milk quickly.

Tips 1

- Only store the breast milk collected by sterilized breast pump.
- Do not freeze the breast milk repeatedly
- Do not mix fresh breast milk with frozen breast milk
- Be sure to discard the remaining breast milk after feeding

Tips 2

- Remember! Always check the milk temperature before feeding.
- Never heat breast milk in microwave as this may destroy valuable nutrients and antibodies.
- Uneven heating may lead to local overheating of the milk, which may cause burns.
- Do not immerse a frozen bottle in boiling water to avoid crack or split.

Note

This product has a built-in lithium battery. When the battery indicator flashes in red, we recommend that you use a certified 5V = 1A adapter to charge the pump motor as soon as possible.